


# VivaStorage™ – Data Management

VivaStorage is a secure way to store and backup all clinical, dermoscopic and confocal data to a professional data center. The data is stored according to highest safety standards and can only be accessed via secure login from any computer connected to the internet.



## ■ Advantages at a Glance

- ▶ No more local data backup required
- ▶ Never run out of storage
- ▶ Long-term archiving of data
- ▶ Automatic data transfer possible
- ▶ Customizable setup to match your needs
- ▶ Access your data from any computer with internet access

## ■ Technical Data – MAVIG Server

- Professional server, hosted by Hartl EDV GmbH & Co KG
- State of the art server network
- Data safety ensured by continuous backup
- High speed access from anywhere
- Modern data center ISO/IEC 27001:2013 certified
- Server located in Germany

*Data Center / Server  
Storage & Exchange*



## ■ Data Safety and Security

- All data is encrypted during transfer using HTTPS protocol
- All services are compliant with EU General Data Protection Regulation (GDPR) privacy and security rules

## ■ Local IT Requirements

- Minimum internet connection: VDSL 25 (upload speed: > 5 Mbit/s)
- Local firewall allowing outgoing HTTPS traffic via port 443

### **MAVIG GmbH** **VivaScope Systems**

Stahlgruberring 5  
81829 Munich, Germany  
Phone +49 (0) 89 / 420 96 - 280  
Email [info@vivascope.eu](mailto:info@vivascope.eu)

[www.vivascope.eu](http://www.vivascope.eu)

